## ENDOCRINOLOGY

# Copeptin aids in the diagnosis of polyuria-polydipsia syndrome



## What is Copeptin?

Copeptin and vasopressin (AVP), often known as antidiuretic hormone (ADH), are derived from the same pre-provasopressin precursor molecule. Serum copeptin concentrations are substantially correlated with AVP levels, and both molecules respond equally to stimuli and changes in blood volume. Unlike AVP, which is released in a pulsatile pattern and swiftly removed from plasma, copeptin levels in plasma tend to be relatively steady. Copeptin is also significantly more stable during storage than AVP, allowing for improved accuracy and easier sample handling.<sup>1</sup>

#### What is the Clinical Utility of Copeptin? Polyuria-polydipsia syndrome: Improved differential diagnosis

We understand the challenges in diagnosing polyuriapolydipsia syndrome and distinguishing between cases of primary polydipsia and diabetes insipidus. Direct AVP measurement would be the preferred diagnostic method, but the molecular properties of AVP make such testing complex and unreliable. A single copeptin assay can now distinguish between central diabetes insipidus and nephrogenic diabetes insipidus. This can reduce the burden of the water deprivation test for the majority of patients while also improving patient management.<sup>1</sup>


#### **Pituitary surgery:** Easy monitoring for vasopressin deficiency

Pituitary tumors, due to their specific location, can cause a variety of hormone deficiencies and water metabolism issues. Copeptin allows for simple monitoring of vasopressin deficiency after pituitary surgery. Low copeptin levels within twelve hours of surgery are a reliable predictor of the development of post-operative diabetic insipidus.<sup>1,2</sup>

### Additional diagnostic potential of copeptin

- Traumatic brain injury—reliable follow-up for hypopituitarism<sup>3</sup>
- Hyponatremia—a promising biomarker for differential diagnosis<sup>4</sup>
- Polycystic kidney disease—better assessment of disease progression and management<sup>5</sup>
- Septic shock—fine-tuning of exogenous vasopressin therapy<sup>6</sup>
- Acute myocardial infarction (AMI)—early rule out measurement in combination with troponin<sup>7</sup>





A diagnostic workflow for the differential diagnosis of polyuria-polydipsia syndrome, modified from Christ-Crain M et al.<sup>1</sup>

Test Name	Test No.	Specimen	Collection	Storage	Patient Preparation
Copeptin	010505	Plasma	Separate plasma from cells and transfer to a plastic transport tube	Room temperature	None required

#### References

1. Christ-Crain M, Fenske W. Copeptin in the diagnosis of vasopressin-dependent disorders of fluid homeostasis. *Nat Rev Endocrinol*. 2016;12(3):168-176. doi:10.1038/nrendo.2015.224. 2. Winzeler B, Zweifel C, Nigro N, et al. Postoperative Copeptin Concentration Predicts Diabetes Insipidus After Pituitary Surgery. *J Clin Endocrinol Metab*. 2015;100(6):2275-2282. doi:10.1210/ jc.2014-4527.

3. Glynn N, Agha A. Which patient requires neuroendocrine assessment following traumatic brain injury, when and how? *Clin Endocrinol (Oxf)*. 2013;78(1):17-20. doi:10.1111/cen.12010. 4. Boursier G, Alméras M, Buthiau D, et al. CT-pro-AVP as a tool for assessment of intravascular volume depletion in severe hyponatremia. *Clin Biochem*. 2015;48(10-11):640-645. doi:10.1016/j. clinbiochem.2015.03.013.

5. Meijer E, Bakker SJ, van der Jagt EJ, et al. Copeptin, a surrogate marker of vasopressin, is associated with disease severity in autosomal dominant polycystic kidney disease. *Clin J Am Soc Nephrol.* 2011;6(2):361-368. doi:10.2215/CJN.04560510.

6. Morgenthaler NG, Müller B, Struck J, Bergmann A, Redl H, Christ-Crain M. Copeptin, a stable peptide of the arginine vasopressin precursor, is elevated in hemorrhagic and septic shock. 2007;28(2):219-226. doi:10.1097/SHK.0b013e318033e5da.

7.Lipinski MJ, Escárcega RO, D'Ascenzo F, et al. A systematic review and collaborative meta-analysis to determine the incremental value of copeptin for rapid rule-out of acute myocardial infarction. Am J Cardiol. 2014;113(9):1581-1591. doi:10.1016/j.amjcard.2014.01.436.



Please contact your local account representative for more information, or visit **Labcorp.com.**