

# BRCA 1/2 Next Generation Sequencing Assay Summary

---

## Analyte

Mutation detection within the *BRCA1* and *BRCA2* genes

## Test System

**Methodology:** Next Generation Sequencing

**Instruments:** Illumina MiSeq®

## Assay Principle

PCR Amplification and Next Generation Sequencing

The entire gene coding region of *BRCA1/BRCA2* genes, as well as all flanking noncoding regions, is analyzed by next generation sequencing. The test uses Illumina's TruSeq® Custom Amplicon approach for target enrichment and MiSeq instrument for sequencing. It is a custom-designed pool comprised of 150 oligonucleotide probe pairs for the binding and amplification of products ranging from 220 bp to 280 bp.

## Typical Applications of Assays

According to the National Comprehensive Cancer Network, testing is indicated if one of the features mentioned below is present in the family: Early-age-onset (age < 50 years) breast cancer, including both invasive and ductal carcinoma in situ (DCIS) breast cancers; two breast primaries or breast and ovarian/fallopian tube/primary peritoneal cancer in a single individual or two or more breast primaries or breast and ovarian/fallopian tube/primary peritoneal cancers in close (first-, second-, and third-degree) relatives(s) from the same side of the family; populations at risk (eg, Ashkenazi Jewish); member of a family with a known *BRCA1* or *BRCA2* mutation; any male breast cancer; ovarian/fallopian tube/primary peritoneal cancer at any age.

## Biological Limitations

Samples must be whole blood

## Preanalytical Considerations

**Volume** 7 mL

**Minimum Volume** 3 mL

**Container** Lavender-top (EDTA) tube or yellow-top (ACD) tube

**Storage Instructions** 2°–8°C (preferred), or room temperature (15°–25°C). Specimens should be shipped at ambient temperature.

**Stability** 48 hours at ambient temperature

**Patient Preparation** None

## Assay Characteristics

Variant detection within the *BRCA1* and *BRCA2* genes is based on Integrated Genetics validation data and CMBP transfer validation data (June 2014).

Parameter	Description	Comments
<b>Accuracy</b>	<ul style="list-style-type: none"> <li>• 72 Clinical and 30 Coriell DNAs used during transfer validation</li> <li>• 781 Het (SNV+indels1-40bp)</li> <li>• 445 Hom (all SNVs)</li> <li>• 71 Clinical and 31 Coriell DNAs used during validation</li> <li>• 1780 Het (SNV+indels1-40bp)</li> <li>• 984 Hom (all SNVs)</li> </ul>	<ul style="list-style-type: none"> <li>• Concordance to Sanger: 100% of mutations identified by NGS</li> <li>• Concordance between sites: 100%</li> </ul>
<b>Intraassay Precision</b>	100% on all replicates of all identified variants on six samples	
<b>Interassay Precision</b>	100% on all replicates of all identified variants on six samples	Second tech and second lot verified
<b>Sensitivity</b>	99.9%	
<b>Specificity</b>	Test is 100% specific for the <i>BRCA1</i> and <i>BRCA2</i> genes	Manufacturer and SNPcheck verified
<b>Minimum Coverage Threshold</b>	30X	All samples except one passed the coverage metric.
<b>Reference Interval</b>	N/A; specific variant reported	
<b>Reportable Range</b>	N/A; variants reported from variants of unknown significance to causative variants	17,769 bp total

## Assay Status

- Variant detection within the *BRCA1* and *BRCA2* genes has been validated at Integrated Genetics, a brand of Esoterix Genetic Laboratories, LLC, a wholly owned subsidiary of LabCorp, Westborough, Mass, and subsequently transferred and verified with additional validation at the Center for Molecular Biology and Pathology, Research Triangle Park, NC.
- The assay is analytically validated and is performed in a CLIA-compliant laboratory at the Center for Molecular Biology and Pathology (Research Triangle Park, NC).
- Validation summary is on file.
- Assay turnaround time is 7 to 12 days.

## Limitations

Sequencing does not detect inversions or other rearrangements.

### References

1. American Society of Clinical Oncology. American Society of Clinical Oncology Policy Statement Update: Genetic testing for cancer susceptibility. *J Clin Oncol*. 2003 Jun 15; 21(12):2397-2406. PubMed 12692171
2. Antoniou A, Pharoah PD, Narod S, et al. Average risks of breast and ovarian cancer associated with *BRCA1* or *BRCA2* mutations detected in case Series unselected for family history: A combined analysis of 22 studies. *Am J Hum Genet*. 2003 May; 72(5):1117-1130. PubMed 12677558
3. Chen S, Parmigiani G. Meta-analysis of *BRCA1* and *BRCA2* penetrance. *J Clin Oncol*. 2007 Apr 10; 25(11):1329-1333. PubMed 17416853
4. Engert S, Wappenschmidt B, Betz B, et al. MLPA screening in the *BRCA1* gene from 1,506 German hereditary breast cancer cases: Novel deletions, frequent involvement of exon 17, and occurrence in single early-onset cases. *Hum Mutat*. 2008 Jul; 29(7):948-958. PubMed 18431737
5. Ferla R, Calò V, Cascio S, et al. Founder mutations in *BRCA1* and *BRCA2* genes. *Ann Oncol*. 2007 Jun; 18(Suppl 6):vi93-vi98. PubMed 17591843
6. Ford D, Easton DF, Stratton M, et al. Genetic heterogeneity and penetrance of the *BRCA1* and *BRCA2* genes in breast cancer families. The Breast Cancer Linkage Consortium. *Am J Hum Genet*. 1998 Mar; 62(3):676-689. PubMed 9497246
7. Herman S, Varga D, Deissler HL, Kreienberg R, Deissler H. Medium-sized deletion in the *BRCA1* gene: Limitations of Sanger sequencing and MLPA analyses. *Genet Mol Biol*. 2012 Jan; 35(1):53-56. PubMed 22481874
8. Janavicius R. Founder *BRCA1/2* mutations in the Europe: Implications for hereditary breast-ovarian cancer prevention and control. *EPMA J*. 2010 Sep; 1(3):397-412. PubMed 23199084
9. Kote-Jarai Z, Leongamornlert D, Saunders E, et al. *BRCA2* is a moderate penetrance gene contributing to young-onset prostate cancer: Implications for genetic testing in prostate cancer patients. *Br J Cancer*. 2011 Oct 11; 105(8):1230-1234. PubMed 21952622
10. National Comprehensive Cancer Network. Clinical practice guidelines in oncology, genetic/familial high-risk assessment: breast and ovarian. Available at: www.nccn.org. 2010. Accessed May 29, 2011
11. Roy R, Chun J, Powell SN. *BRCA1* and *BRCA2*: Different roles in a common pathway of genome protection. *Nat Rev Cancer*. 2011 Dec 23; 12(12):68-78. PubMed 22193408
12. Tai YC, Domchek S, Parmigiani G, Chen S. Breast cancer risk among male *BRCA1* and *BRCA2* mutation carriers. *J Natl Cancer Inst*. 2007 Dec 5; 99(23):1811-1814. PubMed 18042939
13. Thompson D, Easton DF; Breast Cancer Linkage Consortium. Cancer incidence in *BRCA1* mutation carriers. *J Natl Cancer Inst*. 2002 Sep 18; 94(18):1358-1365. PubMed 12237281
14. Tulinius H, Olafsdottir GH, Sigvaldason H, et al. The effect of a single *BRCA2* mutation on cancer in Iceland. *J Med Genet*. 2002 Jul; 39(7):457-462. PubMed 12114473



www.LabCorp.com