



# THYROID TESTING

ASSESSING THYROID DISEASE  
IN YOUR PATIENTS



# THYROID FUNCTION

LabCorp's expertise in thyroid testing provides clinicians with a comprehensive portfolio for their thyroid needs.

## SCREENING FOR THYROID FUNCTION IS IMPORTANT.

It is estimated that 20 million Americans suffer from some form of thyroid disease.<sup>1</sup>

More than 70% of hyperthyroid cases are caused by overproduction of thyroid hormone by the thyroid gland (Graves' disease).<sup>2</sup>

Hypothyroidism accounts for approximately 80% of patients with thyroid disorders.<sup>3</sup>

**Diagnosing whether hyper- or hypothyroidism is caused by an autoimmune disease is critical for patient care and treatment.**

- Graves' disease is the most common form of hyperthyroidism, and diagnostic testing indicators show a presence of thyroid-stimulating immunoglobulin (TSI) and TSH receptor antibodies (TRAb/TBII).<sup>4</sup>
- Hashimoto's thyroiditis is one of the most common forms of hypothyroidism,<sup>5</sup> and it is usually characterized by the development of antithyroid peroxidase (anti-TPO) and/or antithyroglobulin (anti-Tg) autoantibodies.<sup>6</sup>



## LabCorp Offers Tests That Can Assist

in the diagnosis of thyroid disorders, including hyperthyroidism, hypothyroidism, and autoimmune diseases such as Graves' disease and Hashimoto's thyroiditis.



Clients also have direct access to specialized offerings and services through LabCorp's Specialty Testing Group, including Endocrine Sciences and Dianon Pathology. Endocrine Sciences is a research-quality laboratory specializing in highly specific and sensitive endocrine testing services. Dianon Pathology is a leader in providing anatomic pathology services.


## Clinical Application

LabCorp offers several test options to assist with your diagnosis, including:

- Thyroid cascade testing to assist with the diagnosis of thyroid dysfunction.
- A full menu of thyroid tests for diagnosis of thyroid function and autoimmune disease.
- Free thyroxine (T4) and free triiodothyronine (T3) assays by dialysis and HPLC/MS-MS, which provides enhanced accuracy over commonly used analog (automated) methods.<sup>7-10</sup> Dialysis and HPLC/MS-MS may be preferred for evaluating patients who have conditions that impact protein binding capacity, including those who:
  - Have congenital absence of thyroxine-binding protein (TBG).<sup>7</sup>
  - Are pregnant.<sup>11</sup>
  - Are taking oral contraceptives or undergoing hormone therapy.<sup>12</sup>
  - Are taking antipsychotic medications.<sup>13</sup>
  - Have been diagnosed with a malignancy or other critical illness.<sup>14</sup>

## LabCorp Thyroid Cascade

The panel is based on a cascade algorithm that selects specific assays based on the results of previously performed tests, which are necessary to arrive at the most appropriate laboratory diagnosis.



**Note:** Please refer to LabCorp's Thyroid Cascade Technical Review (L964) for additional information regarding interpretation and references for thyroid cascade testing.

# Thyroid Cancer Testing

LabCorp offers comprehensive testing for the diagnosis and monitoring of thyroid cancer.

Detection and diagnosis of thyroid cancer is important.

- Thyroid nodules are common, particularly in older adults. While fewer than 1 in 10 adults have palpable thyroid nodules, when thyroid ultrasonography is performed up to half of adults examined are found to have nodules.<sup>15</sup> About 1 in 20 thyroid nodules is cancerous.<sup>16</sup>
- A thyroid biopsy using fine needle aspiration can differentiate between malignant and benign nodules.
- When diagnosed and treated, common thyroid tumors present an 80% to 90% survival rate at 10 years, with survival rates of greater than 97% for younger patients who are treated appropriately.<sup>16</sup>

LabCorp offers serum calcitonin testing for patients with suspected medullary thyroid carcinoma (MTC)<sup>17</sup>

- Serum calcitonin testing is useful for the detection and confirmation of C-cell hyperplasia (the precursor of MTC) as well as a tumor marker for diagnosis and management of MTC.
- Preoperative serum calcitonin is reported to roughly correlate with tumor weight or extent of disease; therefore, postoperative levels also have prognostic application.

LabCorp also offers molecular blood testing for RET gene mutations.

- Testing includes mutation analysis by sequencing of exons 10, 11, 13, 14, 15, and 16 of the RET proto-oncogene.
- Mutations found in these exons have been associated with patients that develop multiple endocrine neoplasia type 2 (MEN 2) and/or familial medullary thyroid carcinoma (FMTC). Testing for RET germline mutations is recommended in patients with a family or personal history consistent with MEN 2 or FMTC.<sup>17</sup>
- Known mutation testing of RET mutations is also available for patients with family members who have previously identified mutations.

LabCorp and Dianon Pathology, a member of the LabCorp Specialty Testing Group, offer several test options for detecting and diagnosing thyroid cancer using fine needle aspiration (FNA) biopsies.

## Fine Needle Aspirate Collection Kits

**LabCorp and Dianon Pathology both offer fine needle aspirate collection kits that employ a space-saving design to assist with efficient biopsy collection. The collection kits include:**

- 8 slides with fixative containers and 1 CytoLyt® vial.
- Convenient design to hold specimen containers in place during the aspiration procedure.
- Specialized thyroid FNA kits with an option that includes a vial of RNARetain® for molecular testing



### Dianon Test Kit

Dianon Thyroid FNA Kit  
with RNARetain®

If FNA results are indeterminate, LabCorp and Dianon Pathology offer ThyGenX® with reflex to ThyraMIR™. These laboratories also offer an option for the Thyroid FNA test with indeterminate reflex to ThyGenX® only, and utilize the Bethesda system nomenclature for thyroid FNA cytology results.

Dianon Pathology also provides the following cytopathology services, including:

- Specialized endocrine pathology requisition
- Full-color reports with photomicrographs
- Dedicated cytopathology staff with expertise that includes thyroid biopsies

**ThyGenX® and ThyGenX® with reflex to ThyraMIR™ are performed by Interpace Diagnostics as a send-out from LabCorp.**

- ThyGenX® includes markers for BRAF, KRAS, HRAS, NRAS, PIK3CA, RET/PTC1, RET/PTC3, and PAX8/PPARGgamma
- ThyraMIR™ includes 10 miRNA markers and occurs if ThyGenX® is negative or not fully indicative of malignancy
- Combined negative predictive value (likelihood negative result is truly benign) found to be 94%<sup>22</sup>
- Combined positive predictive value (likelihood positive result is malignant) found to be 74%<sup>22</sup>
- Predicted to result in up to 85% reduction in unnecessary diagnostic surgeries<sup>22</sup>

# Thyroid Cancer Monitoring

Once thyroid cancer has been diagnosed and treated, patients must be closely monitored for cancer recurrence.

- Thyroglobulin is a protein secreted only by thyroid tissue.
- After thyroidectomy, thyroglobulin levels are commonly used to detect recurrence of thyroid cancer.<sup>18</sup>
- Patients with MTC are also monitored with calcitonin and CEA testing as recommended by American Thyroid Association guidelines.<sup>17</sup>


LabCorp offers enhanced sensitivity for thyroglobulin and antithyroglobulin testing to monitor for thyroid cancer recurrence.

- Serum thyroglobulin (Tg) is primarily used in the postoperative management of differentiated thyroid cancer (DTC).
- Thyroglobulin antibody (TgAb) is detected in an estimated 25% of patients with DTC.<sup>19</sup> In those patients, there is a risk of interference with Tg measurement using immunometric (IMA) methods that can lead to false-negative (inappropriately low or undetectable) Tg results.<sup>4,19,20</sup> Even low antibody concentrations can interfere with Tg measurements.<sup>4</sup>
- LabCorp's Thyroglobulin With Antithyroglobulin Antibody test offers a dual assay strategy for Tg in an effort to minimize the potential effect of TgAb interference.
  - Specimens are tested for TgAb using a sensitive IMA.
  - Specimens with TgAb below the detectable limit (<1.0 IU/mL) are tested for Tg by sensitive second-generation IMA.
  - Specimens with any measurable TgAb levels ( $\geq 1.0$  IU/mL) are tested for Tg by radioimmunoassay (RIA), which is less prone to interference by TgAb.<sup>4</sup>
- Clinicians who require extended (1 year) specimen storage may order LabCorp's Comprehensive Thyroglobulin Profile, which is performed by Endocrine Sciences, a member of the LabCorp Specialty Testing Group.

LabCorp offers thyroglobulin testing for lymph node aspirate diluted in saline.

- Aspirate material from lymph nodes can be tested for the presence of thyroglobulin if there is suspicion that thyroid cancer has spread to the lymphatic system.<sup>21</sup>
- The lymph node aspirate is collected and washed into a 1 mL saline solution, and the saline solution is tested for the presence of thyroglobulin.

## Thyroglobulin With Antithyroid Antibody Test



## Lymph Node Aspirate Collection Kit

A convenient lymph node aspirate collection kit is available that includes instructions and a vial containing 1 mL of saline solution.



# Thyroid Testing Services

## Thyroid Function: Testing for Hyper- and Hypothyroidism

Test Name	Number	Methodology	Specimen	Container	Storage
<b>Thyroid Cascade Profile: TSH with automatic reflex (as diagnostically warranted) to FT<sub>4</sub>, FT<sub>3</sub>, and/or TPO antibodies</b>	<b>330015</b>	ECLIA	2.0 mL serum Minimum: 1.0 mL	Red-top tube or gel-barrier tube	Refrigerate
<b>Thyroxine (T<sub>4</sub>)</b>	<b>001149</b>	CEDIA	Serum: 1 mL (adult), 0.8 mL (pediatric) Minimum: 0.5 mL (adult), 0.3 mL (pediatric)	Red-top tube or gel-barrier tube	Refrigerate
<b>Thyroid-stimulating Hormone (TSH)</b>	<b>004259</b>	ECLIA	0.8 mL serum Minimum: 0.3 mL	Red-top tube or gel-barrier tube	Refrigerate
<b>Thyroxine (T<sub>4</sub>) Free, Dialysis/Mass Spectrometry<sup>++</sup></b>	<b>501902*</b>	Direct dialysis mass spectrometry; HPLC/MS	1.0 mL serum <b>or</b> plasma Minimum: 0.5 mL	Red-top tube or lavender-top (EDTA) tube	<b>Freeze</b>
<b>Thyroxine (T<sub>4</sub>), Free, Direct, Serum</b>	<b>001974</b>	ECLIA	0.8 mL serum Minimum: 0.3 mL	Red-top tube or gel-barrier tube	Refrigerate
<b>Thyroxine-binding Globulin (TBG), Serum</b>	<b>001735</b>	ICMA	0.5 mL serum Minimum: 0.3 mL	Red-top tube or gel-barrier tube	Refrigerate
<b>Triiodothyronine (T<sub>3</sub>)</b>	<b>002188</b>	ECLIA	0.8 mL serum Minimum: 0.3 mL	Red-top or gel-barrier tube	Refrigerate
<b>Triiodothyronine (T<sub>3</sub>), Free, Serum</b>	<b>010389</b>	ECLIA	0.8 mL serum Minimum: 0.3 mL	Red-top tube or gel-barrier tube	Refrigerate
<b>Triiodothyronine, Free (FT<sub>3</sub>), Dialysis and LC-MS/MS<sup>++</sup></b>	<b>503600*</b>	Equilibrium dialysis and HPLC/MS-MS	1.0 mL serum Minimum: 0.3 mL	Red-top tube, gel-barrier tube, lavender-top (EDTA) tube, or green-top (sodium heparin) tube	<b>Freeze</b> (preferred) or refrigerate

## General Autoimmune Screen: Testing for Hyper- and Hypothyroidism

<b>Thyroglobulin Antibody</b>	<b>006685</b>	ICMA	1.0 mL serum	Red-top tube or gel-barrier tube	Room Temperature
<b>Thyroid Antibodies</b>	<b>006684</b>	See individual test descriptions	1.0 mL serum	Red-top tube or gel-barrier tube	Refrigerate
<b>Thyroid Peroxidase (TPO) Antibodies</b>	<b>006676</b>	ECLIA	0.8 mL serum Minimum: 0.3 mL	Red-top tube or gel-barrier tube	Refrigerate

## Graves' Disease Autoimmune Screen

<b>Thyroid-Stimulating Immunoglobulin (TSI)</b>	<b>140749</b>	Chinese hamster ovary cell line transfected with thyrotropin receptor and a luciferase reporter gene	3.0 mL serum Minimum: 0.3 mL	Red-top tube or gel-barrier tube	Refrigerate
<b>TSH Receptor Antibody (TRAb/TBII)</b>	<b>500538*</b>	Binding inhibition assay	1.0 mL serum Minimum: 0.3 mL	Red-top tube or gel-barrier tube	Ambient (same day) or <b>freeze</b>

Visit the online Test Menu at [www.LabCorp.com](http://www.LabCorp.com) for full test information, including CPT codes and current specimen collection requirements.

\* Testing performed at Endocrine Sciences.

++ Free T<sub>4</sub> and free T<sub>3</sub> by dialysis and HPLC/MS-MS should be used for patients known to have abnormal binding proteins due to pregnancy, hormone replacement, or critical illnesses.



## Thyroid Cancer Screen: Testing for Diagnosis of Thyroid Cancer

Test Name	Number	Methodology	Specimen	Container	Storage
Calcitonin (Thyrocalcitonin)	004895	ICMA	1.0 mL serum Minimum: 0.4 mL	Red-top tube or gel-barrier tube	Freeze
Fine Needle Aspiration Cytology	009001	Morphologic analysis	Aspirated material Recommend using LabCorp Collection kit: Catalog N° FNAK10	Slide(s); Coplin jar(s)	Refrigerate
MEN2: RET Gene Sequencing (for hereditary thyroid cancer)	504008*	PCR, sequencing	4.0 mL whole blood Minimum: 1.0 mL	Lavender-top (EDTA) tube	Ambient or refrigerate

## Thyroid Cancer Monitoring: Testing for Recurrence of Thyroid Cancer

Calcitonin (Thyrocalcitonin)	004895	ICMA	1.0 mL serum Minimum: 0.4 mL	Red-top tube or gel-barrier tube	Freeze
Carcinoembryonic Antigen (CEA)	002139	ECLIA	0.8 mL serum Minimum: 0.3 mL	Red-top tube or gel-barrier tube	Refrigerate
Thyroglobulin Antibody and Thyroglobulin, IMA or LC/MS-MS	042045	IMA or LC/MS-MS	3 mL serum Minimum: 2 mL	Red-top tube or gel-barrier tube	Room temperature
Thyroglobulin Antibody and Thyroglobulin, IMA or RIA	042060	IMA or RIA	3 mL serum Minimum: 2 mL	Red-top tube or gel-barrier tube	Room temperature
Thyroglobulin, Lymph Node Aspirate	502380*	ICMA	Lymph node aspirate in 1.0 mL saline Recommend using LabCorp Collection kit: Catalog N° 38621G	Lymph Node Collection Kit and saline vial	Freeze

Visit the online test menu at [www.LabCorp.com](http://www.LabCorp.com) for full test information, including CPT codes and current specimen collection requirements.

\* Testing performed at Endocrine Sciences.

++ Free T<sub>4</sub> and free T<sub>3</sub> by dialysis and HPLC/MS-MS should be used for patients known to have abnormal binding proteins due to pregnancy, hormone replacement, or critical illnesses.

## References

1. American Association for Clinical Chemistry. Thyroid diseases. [Lab Tests Online Web site]. December 20, 2009. Available at: <http://www.labtestsonline.org/understanding/conditions/thyroid.html>. Accessed February 8, 2010.
2. American Thyroid Association. Hyperthyroidism. Available at: [http://www.thyroid.org/patients/patient\\_brochures/hyperthyroidism.html](http://www.thyroid.org/patients/patient_brochures/hyperthyroidism.html). Accessed December 22, 2009.
3. Hollowell JG, Staehling NW, Flanders WD, et al. Serum TSH, T4, and thyroid antibodies in the United States population (1988 to 1994): National health and nutrition examination survey (NHANES III). *J Clin Endocrinol Metab.* 2002;87:489-499.
4. Demers LM, Spencer CA. *Laboratory Support for the Diagnosis of Thyroid Disease*. Washington, DC: The National Academy of Clinical Biochemistry; 2002. Volume 13.
5. American Thyroid Association. Hypothyroidism. Available at: [http://www.thyroid.org/patients/patient\\_brochures/hypothyroidism.html](http://www.thyroid.org/patients/patient_brochures/hypothyroidism.html). Accessed December 22, 2009.
6. American Association for Clinical Chemistry. Thyroid antibodies. [Lab Tests Online Web site]. December 20, 2009. Available at: [http://labtestsonline.org/understanding/analytes/thyroid\\_antibodies/test.html](http://labtestsonline.org/understanding/analytes/thyroid_antibodies/test.html). Accessed December 22, 2009.
7. Yue B, Rockwood AL, Sandrock T, La'ulu SL, Kushnir MM, Meikle AW. Free thyroid hormones in serum by equilibrium dialysis and online solid-phase extraction—liquid chromatography/tandem mass spectrometry. *Clin Chem.* 2008;54:642-651.
8. Soldin SJ, Soukhova N, Janicic N, Jonklaas J, Soldin OP. The measurement of free thyroxine by isotope dilution tandem mass spectrometry. *Clinica Chimica Acta.* 200; 358:113-118.
9. Nelson JC, Weiss RM, Wilcox RB. Underestimates of serum free thyroxine (T4) concentrations by free T4 immunoassays. *J Clin Endocrinol Metab.* 1994;79:76-79.
10. Nelson JC, Tomel RT. Direct determination of free thyroxine in undiluted serum by equilibrium dialysis/radioimmunoassay. *Clin Chem.* 1988;34:1737-1744.
11. Kahric-Janjic N, Soldin SJ, Soldin OP, West T, Gu J, Jonklaas J. Tandem mass spectrometry improves the accuracy of free thyroxine measurements during pregnancy. *Thyroid.* 2007;17:303-311.
12. Mazer N. Interaction of estrogen therapy and thyroid hormone replacement in postmenopausal women. *Thyroid.* 2004;14:S27-34.
13. Surks MI, DeFesi CR. Normal serum free thyroid hormone concentrations in patients treated with phenytoin or carbamazepine. A paradox resolved. *JAMA.* 1996; 275:1495-1498.
14. De Groot L. Dangerous dogmas in medicine: the nonthyroidal illness syndrome. *J Clin Endocrinol Metab.* 1999;84:151-164.
15. Frates MC, Benson CB, Charboneau JW, et al. Management of thyroid nodules detected at US: Society of Radiologists in Ultrasound consensus conference statement. *Radiology.* 2005;237:794-800.
16. American Cancer Society. What is thyroid cancer? Available at: [http://www.cancer.org/docroot/CRI/content/CRI\\_2\\_4\\_1X\\_What\\_is\\_thyroid\\_cancer\\_43.asp?sitearea=](http://www.cancer.org/docroot/CRI/content/CRI_2_4_1X_What_is_thyroid_cancer_43.asp?sitearea=). Accessed February 9, 2010.
17. Kloos RT, Eng C, Evans DB. Medullary thyroid cancer: Management guidelines of the American Thyroid Association. *Thyroid.* 2009;19:565-613.
18. Spencer CA, Bergoglio LM, Kazarosyan M, Fatemi S, LoPresti JS. Clinical impact of thyroglobulin (Tg) and Tg autoantibody method differences on the management of patients with differentiated thyroid carcinomas. *J Clin Endocrinol Metab.* 2005 Oct; 90(10):5566-5575.
19. Cooper DS, Doherty GM, Haugen BR, Kloos RT, et al. Revised American Thyroid Association management guidelines for patients with thyroid nodules and differentiated thyroid cancer. *Thyroid.* 2009;19:1-48.
20. Gharib H, Papini E, Paschke R, et al. American Association of Clinical Endocrinologists, *Associazione Medici Endocrinologi*, and European Thyroid Association medical guidelines for clinical practice for the diagnosis and management of thyroid nodules. *Endocr Pract.* 2010;16(3):468-475.
21. Borel AL, Boizel R, Faure P, et al. Significance of low levels of thyroglobulin in fine needle aspirates from cervical lymph nodes of patients with a history of differentiated thyroid cancer. *Euro J Endocrinol.* 2008;158:691-698.
22. Labourier E, Shifrin A, Busseniers AE, Lupu MA, Manganelli ML, Andrus B, Wylie D, Beaudenon-Huibregtse S. Molecular testing for miRNA, mRNA and DNA on fine needle aspiration improves the preoperative diagnosis of thyroid nodules with indeterminate cytology. *J Clin Endocrinol Metab.* 2015;100(7):2743-2750. PubMed 25965083

## Contact Us

Please contact your local account representative for more information, or visit [www.LabCorp.com](http://www.LabCorp.com).



[www.LabCorp.com](http://www.LabCorp.com)