# **GenoSure Archive**°

Cell-associated viral DNA assay for HIV-1 protease, reverse transcriptase, and integrase inhibitors

### Introduction

Human immunodeficiency virus type 1 (HIV-1) was identified in 1983 as the cause of acquired immune deficiency syndrome (AIDS). During the decade between 1985 and 1995, more than 100,000 people in the United States lost their lives to HIV.<sup>1</sup> In March of 1987, the US Food and Drug Administration (FDA) approved the first antiretroviral (ARV) drug azidothymidine (AZT; zidovudine) for the treatment of HIV/AIDS, and more than 30 ARV drugs comprising six drug classes have been approved to date.<sup>2</sup>

By design, most ARV drugs target the activity of the viral enzymes required by HIV-1 to complete its lifecycle: nucleoside reverse transcriptase inhibitors (NRTIs) and nonnucleoside reverse transcriptase inhibitors (NNRTIs) target the activity of the reverse transcriptase enzyme; integrase inhibitors (INIs) target the activity of the integrase enzyme; and protease inhibitors (PIs) target the activity of the protease enzyme. ARV drugs, however, do not target, disrupt, or delete integrated copies of viral DNA (provirus) that are archived within the host cell genome during replication. Thus HIV remains an incurable, yet manageable, disease.


The plurality of approved ARV drugs combined with therapeutic advances, such as improved efficacy, lower pill burden, less frequent dosing, and reduced side effects, have resulted in more patients with complete and sustained suppression of virus replication. This has redefined HIV/AIDS from a terminal illness to a chronic, manageable disease.<sup>3</sup> Advancements in ARV therapy have increased life expectancy for HIV-1-infected patients in the US and Canada from 36 years in 2000 to 51 years in 2006.<sup>4</sup> The success of ARV therapy for each patient is predicated on the identification of a combination of ARV drugs that is active against HIV-1. During replication, the virus may select for mutations that confer reduced susceptibility/resistance to the current ARV regimen. The ability of the virus to develop ARV resistance was first documented and described in the early 1990s when the use of AZT as monotherapy was common.<sup>5-7</sup> Drug resistance has since been described for all ARV drugs approved for HIV-1 treatment.<sup>8</sup> Consequently, the Department of Health and Human Services recommends screening for drug resistant virus prior to treatment initiation and prior to retreatment after ARV failure.<sup>9</sup>

Resistance tests are routinely used to identify effective ARV drug options so that suppressive regimens can be administered to each individual patient. ARV drug regimens are typically switched or adjusted when the existing/current regimen fails to maintain complete viral suppression. Increasingly common today is the need for "fine tuning" regimens while the patient's virus remains suppressed. Reasons for these treatment adjustments most often include side effects, adverse events, regimen intolerance, and drug-drug interactions.

Standard resistance testing to guide ARV drug selection cannot be employed in the setting of viral suppression because the currently available assays sample plasma virus and require viral loads ≥500 copies/mL. Consequently, new diagnostic tools are needed to give treatment providers ARV drug susceptibility information in the setting of complete viral suppression (undetectable viral load) or low level viremia (<500 copies/mL). These new diagnostic tools exploit the archived HIV-1 proviral DNA that is embedded in host cells during virus replication. Access to this genetic information is enabled by the lysis of HIV-1-infected cells in whole blood samples. The first of these diagnostic tools offered by Monogram Biosciences<sup>®</sup>, Trofile<sup>®</sup> DNA, was launched in 2010. Trofile® DNA evaluates the envelope (env) region of the HIV-1 genome and is used to determine tropism in the setting of complete virologic suppression to guide CCR5 antagonist administration.<sup>10</sup> An additional tool is now available: GenoSure Archive®, which interrogates the HIV-1 polymerase (pol) region to identify archived mutations. This new assay provides susceptibility assessments for PIs, NRTIs, NNRTIs, and INIs, which enable providers to obtain comprehensive ARV drug susceptibility information in the setting of virologic suppression and low-level viremia. This assay facilitates adjustments for a majority of ARV drug regimens.

### Laboratory Method

**GenoSure Archive** is performed by amplifying cellassociated HIV-1 DNA from infected cells in whole-blood samples and then analyzing the pol region using nextgeneration sequencing (NGS) methods. The assay is analytically validated to identify resistance-associated mutations in the HIV-1 pol gene. **GenoSure Archive** evaluates the full-length protease and integrase coding regions and amino acids 1 to 400 of the HIV-1 reverse transcriptase region. (Figure 1).



Briefly, the **GenoSure Archive** requires one 4 mL tube of whole blood. Total DNA is extracted from whole blood, and nested PCR is performed to amplify specifically the HIV-1 region of interest. The resulting amplification products are purified and prepared for nucleic acid sequencing using next generation sequencing methods. After sequencing, the data are analyzed for the presence of resistance-associated mutations using Monogram's proprietary bioinformatics pipeline. Susceptibility calls of "resistant," "sensitive," or "resistance possible" are assigned. These assignments employ the same genotypic interpretive database used in Monogram's GenoSure PRIme® assay, relying on >100,000 paired genotypic and phenotypic susceptibility determinations.

## **Clinical Application**

Studies evaluating the concordance between the HIV-1 plasma RNA and cellular DNA compartments have demonstrated variability. This discordance is likely due to differences that exist between the viral population currently circulating in the plasma and the proviral DNA archived in infected cells. Of note, ARV drugresistant HIV-1 variants in the context of emerging failure are identified earlier in the plasma relative to peripheral blood mononuclear cells (PBMCs).<sup>11</sup> In contrast, resistant variants may persist longer in the latter compartment,<sup>12</sup> especially in the absence of ongoing drug pressure.

A recent study demonstrated favorable concordance, however, between Monogram's HIV-1 DNA genotyping assay and historical plasma genotypes among enrolled patients.<sup>13</sup> In this study, patients who were fully suppressed on their ARV drug regimens, consisting of a boosted protease inhibitor and two NRTIs, were switched to a single-tablet regimen consisting of rilpivirine, tenofovir, and emtricitabine. The authors concluded that cellular DNA genotyping may have clinical utility in detecting archived mutations in aviremic patients prior to switching ARV regimens.

**GenoSure Archive** is designed to provide HIV-1 ARV drug resistance data when standard resistance testing cannot be performed due to inadequate plasma viral load. The assay provides a useful tool in suppression management, specifically to aid in regimen switches or simplification.<sup>14</sup> Often, ARV regimen changes are needed to address side effects, adverse events, and concerns regarding long-term toxicities. There also may be a desire to streamline an overly complex regimen, implemented during an earlier ARV treatment era, when less ideal options existed. Additionally, in the absence of historical drug resistance data for patients maintaining viral suppression, cell-associated viral DNA resistance testing may provide useful information.<sup>15</sup>

The descriptive information above also applies to HIV-1 DNA Sequencing Protease – Reverse Transcriptase and HIV-1 DNA Sequencing Integrase. HIV-1 DNA Sequencing Protease – Reverse Transcriptase provides susceptibility data for nucleosides, nucleotides, nonnucleosides, and protease inhibitors. HIV-1 DNA Sequencing Integrase provides susceptibility data for integrase inhibitors.

#### **Relevant Assays\***

Test Name	Test N°
Human Immunodeficiency Virus 1 (HIV-1) GenoSure Archive®	551776
Human Immunodeficiency Virus 1 (HIV-1) DNA Sequencing Protease – Reverse Transcriptase	551730
Human Immunodeficiency Virus 1 (HIV-1) DNA Sequencing Integrase	551965
Human Immunodeficiency Virus 1 (HIV-1) GenoSure Archive® Plus Trofile® DNA	552020

\*For the most current information regarding test options, including specimen requirements and CPT codes, please consult the online Test Menu at www.LabCorp.com.

#### References

1. Centers for Disease Control and Prevention. *HIV Surveillance Report* — *Epidemiology of HIV Infection*, 2011; vol 23. Available at: http://www.cdc. gov/hiv/pdf/statistics\_epidemiology\_of\_infection\_through\_2011.pdf. February 2013.

2. US Food and Drug Administration. *Antiretroviral drugs used in the treatment of HIV infection*. Available at: http://www.fda.gov/ForConsumers/byAudience/ForPatientAdvocates/HIVandAIDSActivities/ucm118915.htm. Accessed July 24, 2014.

3. Deeks S, Lewin S, Havlir D. The end of AIDS: HIV infection as a chronic disease. *Lancet*. 2012 Nov 2; 382(9903):1525–1533.

4. Samji H, Cescon A, Hogg RS, et al; for The North American AIDS Cohort Collaboration on Research and Design (NA-ACCORD). Closing the gap: Increases in life expectancy among treated HIV-positive individuals in the United States and Canada. *PLoS One.* 2013 Dec 18; 8(12): e81355.

5. Larder BA, Darby G, Richman DD. HIV with reduced sensitivity to zidovudine (AZT) isolated during prolonged therapy. *Science*. 1989 Mar 31; 243:1731–1734.

6. Wainberg MA, Rooke R, Tremblay M, et al. Clinical significance and characterization of AZT-resistant strains of HIV-1. *Can J Infect Dis*. 1991 Spring; 2(1): 5–11.

7. Larder BA, Kemp SD. Multiple mutations in HIV-1 reverse transcriptase confer high-level resistance to zidovudine (AZT). *Science*. 1989 Dec 1; 246(4934):1155-1158.

8. Wensing AM, Calvez V, Günthard HF, et al. 2014 update of the drug resistance mutations in HIV-1. *Top Antivir Med.* 2014 Jun-Jul; 22(3): 642-650.

9. Panel on Antiretroviral Guidelines for Adults and Adolescents. Guidelines for the use of antiretroviral agents in HIV-1-infected adults and adolescents. Department

of Health and Human Services. Available at: http://aidsinfo.nih.gov/contentfiles/ lvguidelines/AdultandAdolescentGL.pdf. Accessed July 23, 2014.

10. Toma J, Frantzell A, Cook J, et al. Phenotypic determination of HIV-1 coreceptor tropism using cell-associated DNA derived from blood samples. Poster presented at: Interscience Conference on Antimicrobial Agents and Chemotherapy; September 11-15, 2010; Boston. Poster H-921.

11. Simmonds P, Zhang LQ, McOmish F, Balfe P, Ludlam CA, Brown AJ. Discontinuous sequence change of human immunodeficiency virus (HIV) type 1 env sequences in plasma viral and lymphocyte-associated proviral populations in vivo: Implications for models of HIV pathogenesis. *J Virol*. 1991 Nov; 65(11):6266-6276.

12. Vandamme AM, Camacho RJ, Ceccherini-Silberstein F, et al. European recommendations for the clinical use of HIV drug resistance testing: 2011 update. *AIDS Rev.* 2011 Apr-Jun; 13(2):77-108.

13. White KL, Toma J, Napolitano LA, et al. Genotypic analyses of pre-existing HIV-1 drug resistance in proviral HIV-1 DNA from PBMCs in suppressed patients switching to RPV/FTC/TDF. Poster presented at: International Workshop on HIV & Hepatitis Virus Drug Resistance and Curative Strategies; June 4-8, 2013; Toronto, Canada. Poster 57. (http://www.natap.org/2013/ResisWksp/ResisWksp\_12.htm).

14. Vandekerckhove LPR, Wensing AMJ, Kaiser R, et al. European guidelines on the clinical management of HIV-1 tropism testing. *Lancet Infect Dis.* 2011 May; 11(5): 394-407.

15. Booth CL, McCormick A, García-Díaz A, et al. Feasibility of testing and detection of HIV-1 drug resistance in proviral DNA. Presented at: The 7th Romanian National HIV/AIDS Congress and The 2nd Central European HIV Forum; 29-31 May 2014; Sibiu, Romania. *BMC Infectious Diseases*. 2014; 14(Suppl 4):O25. GENOSURE ARCHIVE | A TECHNICAL REVIEW



