

Labcorp *AccuDraw*®

Combining on-screen visual cues with step-by-step tutorial, Labcorp *AccuDraw* is an innovative tool designed to support phlebotomists throughout the specimen collection process.

AccuDraw provides a detailed set of specimen collection and handling instructions for your office staff, helping reduce the need to call the lab for collection requirements and allowing staff to spend more time on patient care.

Using *AccuDraw* is simple and straightforward. Prior to the collection, your staff enters the test number(s) to be collected on the labcorp.com *AccuDraw* page. In moments, advanced collection instructions with color-matched, appropriately scaled images of the collection tubes needed will appear on screen. These instructions are customized to the specific physician order and linked to the laboratory's requirements to help ensure accuracy and precision when collecting specimens in your office.

Using *AccuDraw*:

Benefits the physician office and phlebotomist

- **Precision and accuracy**
Visual cues are customized to the specific physician test order and linked to the laboratory's current test requirements
- **High-quality specimens**
Alignment between specimen collection, laboratory requirements and your office helps to ensure the highest quality specimen is sent to the laboratory for testing
- **Time efficiencies**
Reduce the amount of time spent speaking to the lab about specimen requirements
- **Reduction in repeat draws**
AccuDraw's collection, handling and shipping guidance has been proven to reduce significantly the need for patient redraws and repeat visits

Benefits the laboratory

- **Increased efficiencies**
Reduction in need to contact physician about specimen collected or the need for a patient redraw
- **Improved specimen integrity**
Following collection instructions per *AccuDraw* instructions ensures high-quality specimen is received in the lab

Benefits the patients

- Minimize repeat visits for redraws
- Highest quality specimen submitted to lab for testing
- Collection of optimal number of specimens required for testing

Access Labcorp *AccuDraw*

Visit labcorp.com and access Labcorp Link™. You will need to register for access to the Labcorp Link e-services portal.

An *AccuDraw* web plug-in is also available to support integration within existing practice management systems. For additional information on the *AccuDraw* web plug-in, please contact your local Labcorp IT contact.



On-screen visual cues are combined with detailed collection instructions to help ensure specimen collection accuracy within your office.

