



LAB-IN-A-BOX®

A lab solution for home healthcare providers



Kit includes:




Original/
outer box



Inner Container




Test request form (TRF)
with bar code labels



Dual temp packs



Biohazard bag with
tubes



Luer adapter,
vacutainer holder



OSHA compliant
retractable straight
needle



FedEx Pak*

*FedEx Pak is contained only in ordered item “HH5000” and used when shipping Lab-in-a-Box via FedEx

Using Lab-in-a-Box® - Easy as “1, 2, 3”

Lab-in-a-Box provides all of the supplies required for easy blood collection, in 3 simple steps:



Step 1: Draw it!

Perform collection



Step 2: Pack it!

Pack specimens
into Lab-in-a-Box
with temp packs per
instructions, along with
test request form



Step 3: Drop it!

Choose a convenient
delivery option




Features

- Lab-in-a-Box is a proprietary kit designed to serve the needs of home healthcare providers
- Simplifies the blood collection process
- Requires no spinning before drop-off
- Maintains high specimen quality

Laboratory values that can help evaluate the condition of the home healthcare patient

Delivery options - convenient & simple

Delivery and specimen stability




Optimal

Option 1: **Labcorp Patient Service Center (PSC)**

- Convenient and easy
- More than 2000 PSCs nationwide
- Rapid turnaround time

Optimal specimen stability

For specimens that are delivered to a PSC during regular business hours, Labcorp will typically report results for Lab-in-a-Box by the next morning.



Ideal

Option 2: **Labcorp Drop Boxes**

- Convenient and easy
- Strategically placed
- Rapid turnaround time

Ideal specimen stability

This option allows Labcorp to receive the specimen sooner and typically report results for Lab-in-a-Box the next morning.

Note: This offering requires additional logistical planning prior to use.



Acceptable

Option 3: **Federal Express (FedEx®)**

- Convenient
- Only available for certain tests
- When turnaround time is flexible

Acceptable specimen stability for certain tests

This option is only available for certain tests.

Easy-to-read, electronic results

Labcorp has the flexibility to mix and match connectivity solutions to provide efficient and effective communication between Labcorp and the home healthcare provider. Results are available through the following communication channels:

Labcorp Beacon Link® online results

- View, search and trend results
- Easily share results (print, fax, e-mail)
- View trends
- Centralize important notes
- Automated print
- Mobile access to lab results

Fax

- Ability to send fax to multiple providers
- Resolve patient issues

Benefits for all members of your home healthcare team



Cost/time savings

- Lab supplies provided
- Multiple delivery options for added convenience
- Access to online results
- Helps reduce lab redraws



Results

- Timeliness
- Easy to access (fax/online)
- Trending capability



Quality

- The same tests used by physicians daily in their offices and clinics



Nationwide locations

- Standardized operations and processes
- More than 2000 PSCs



Managed care

- Contracts nationwide



Service excellence

- Superior laboratory services

To set up an account

1. Email homehealthcare@labcorp.com with the title "Request a new account"
2. You will receive forms and directions to establish an account
3. Complete all the forms with the requested information, sign and return

Labcorp will provide you with your new account number and information.